

## Product Change Notification

Product Obsolete: DataPro 2.5  
New Cross Reference: DataPro 2.5E  
Effective Date: November 30, 2015

### Description:

DataPro 2.5 is optional facilities monitoring software used with portable Climet particle counters. Data security is compliant with FDA 21 CFR Part 11; and the software provides pass/fail conformance reports to ISO 14644-1:1999, EU GMP Annex 1, and FS 209E standards. The software also provides a scheduling feature, along with graphs and trends.

Effective per the above date, DataPro 2.5 is being obsoleted and replaced by DataPro 2.5E.

All support for the obsoleted DataPro 2.5 will end immediately on all versions purchased over 1 year ago. Otherwise, support will terminate within 12 months of original purchase date.

### Reason:

DataPro 2.5E offers much the same functionality, PLUS added functionality including the following:

- Ethernet and WiFi sample data transfer
- User uploadable floor map
- Ability to monitor multiple CI-x5x portable particle counters on a network.
- Reports include the *new ISO 14644-1:2015 standard*, GMP reports, and Fed Standard 209e

### Impact on Current Customers:

If a customer owns discontinued Climet equipment, the new DataPro 2.5E software may or may not work with discontinued instrumentation. As always, Climet will use best efforts to continue to repair and calibrate obsolete equipment. However, no software support is provided for obsolete instrumentation.

If customer upgrades their operating system or otherwise changes their current computer configuration, the obsoleted DataPro 2.5 software may or may not be compatible with these configuration changes. No support is provided to obsolete software after 12 months of purchase.

**For a limited time, all customers who own DataPro 2.5 may purchase an upgrade of DataPro 2.5E for 50% off the retail price. This offer will expire midnight on 31 December, 2016. Climet may without advance notification, and within its sole and exclusive discretion modify this offer at anytime.**

### **Impact on New Customers:**

No price increase when compared to the older version of the software.

New customers will have the latest software solution including Ethernet and WiFi connectivity, and the ability to generate ISO reports compliant to the 2015 standard.

No impact expected when used with Climet portable particle counters still in manufacture. These models include:

<b>1 CFM Models</b>	<b>50 LPM Models</b>	<b>75 LPM Models</b>	<b>100 LPM Models</b>
CI-150t	CI-450t	CI-750t	
CI-152x	CI-452x	CI-752x	CI-1052x
CI-153	CI-453	CI-753	CI-1053
CI-154	CI-454	CI-754	CI-1054
CI-155	CI-455	CI-755	CI-1055
CI-156	CI-456	CI-756	CI-1056

Climet Application Engineers are committed to assisting Climet customers in making a seamless transition. Please, contact our Sales Team for further assistance.